

5. Roberts, D. L., F. J. Romay, V. A. Marple, "Nozzle Examination Methods for the Next Generation Pharmaceutical Impactor," presented at Drug Delivery to the Lungs XII, The Aerosol Society, London, UK, December 13, 2001.

### **User Experience**

1. Leung, K., E. Louca, M. Gray, G. Tipples, A. L. Coates, "Use of the Next Generation Pharmaceutical Impactor for Particle Size Distribution Measurements of Live Viral Aerosol Vaccines," J. Aerosol Med., 18, 414-26 (2005).
2. Williams, R., A. Sleep, T. Morris, "Next Generation Impactor: Practical Experience of Development and Validation of Metered-Dose Inhaler Methods," Drug Delivery to the Lung XV, The Aerosol Society, London, UK, December 9-10, 2004, pg. 166-69.
3. Beighton, J., A. Majid, L. Wojnarowicz, "The Determination of Aerodynamic Particle Size Distribution from a Metered-Dose Inhaler Using the Next Generation Impactor," Drug Delivery to the Lung XV, The Aerosol Society, London, UK, December 9-10, 2004, pg. 252-55.
4. Mitchell, J. P., "The Next Generation Impactor (NGI): Results from the Evaluation of Prototype Instruments with pMDI-Based Formulations," presented at Drug Delivery to the Lung XI, London, England, December 11-12, 2000.
5. Shrubbs, I. T., "The Next Generation Impactor (NGI): Results from the Evaluation of Prototype Instruments with DPI-Based Formulations," presented at Drug Delivery to the Lung XI, London, England, December 11-12, 2000.

### **Comparisons to Other Sizing Methods**

1. Myrdal, P. B., E. Mogalian, J. P. Mitchell, M. Nagel, C. Wright, B. Kiser, M. Prell, M. Woessner, S. W. Stein, "Application of Heated Inlet Extensions to the TSI 3306/3321 System: Comparison with the Andersen Cascade Impactor and Next Generation Impactor," J. Aerosol Med., 19, 533-42 (2006).

2. Kusmartseva, O., P. R. Smith, D. Morton, "Comparing In-Line Optical and NGI Measurements on a Drug Cloud," Drug Delivery to the Lung XV, The Aerosol Society, London, UK, December 9-10, 2004, pg. 183-86.
3. Kamiya, A., M. Sakagami, M. Hindle, P. R. Byron, "Aerodynamic Sizing of Metered Dose Inhalers: An Evaluation of the Andersen and Next Generation Pharmaceutical Impactors and Their USP Methods," J. Pharm. Sci., 93(7), 1828-37 (2004).
4. Mitchell, J. P., M. W. Nagel, K. J. Wiersema, C. C. Doyle, "Aerodynamic Particle Size Analysis of Aerosols from Pressurized Metered-Dose Inhalers," Pharm. Sci. Tech., 4, article 54, 2003.

#### **Simple or Semi-Automatic Tools for Drug Recovery from the NGI**

1. Roberts, W. L., P. W. Waters, H. Mohammed, R. Williams, B. Gillet, "A Semi-Automated Sample Recovery System for the Next Generation Pharmaceutical Impactor," Respiratory Drug Delivery 2006, Boca Raton, Florida, April 23-27, 2006.
2. Gibbison, F., R. Manley, A. D. Cooper, "Evaluation of the Gentle Rocker as a Laboratory Aid for Fine Particle Dose Determination," Academy of Pharmaceutical Sciences, Bath, UK, March 27-28, 2007.
3. Fransson, K., M. Persson, M. Svensson, "Sample Preparation Tools for the Next Generation Pharmaceutical Impactor," Drug Delivery to the Lung XII, Aerosol Society, London, England, December 13-14, 2001.


## **Appendix A**

### **Major Parts for the Next Generation Pharmaceutical Impactor**

This appendix shows each of the parts of the impactor.




**Figure A-1 – Top Side of NGI Lid**



**Figure A-2 – Bottom Side of NGI Lid Showing the Interstage Passageways**

NGI-0170-6001, Revision B



**Figure A-3 – Top Side of NGI Seal Body with Nozzles and O-Rings**




**Figure A-4 – The Eight Nozzle Pieces of the NGI**

NGI-0170-6001, Revision B




**Figure A-5 – Interference Fit of Nozzle Pieces into Seal Body**




**Figure A-6 – One Large Cup (L) and One Standard Cup (R) – Six Standard Cups and Two Large Cups Comprise a Set**

NGI-0170-6001, Revision B




**Figure A-7 – The Cup Tray**




**Figure A-8 – The Bottom Frame with Handle Attached**



NGI-0170-6001, Revision B




**Figure A-9 – One Set of Cups in a Cup Tray on the Bottom Frame**



**Figure A-10 – Seal Body and Cup Tray Resting on Bottom Frame**

NGI-0170-6001, Revision B



**Figure A-11 – Fully Assembled Impactor in the Open Position**



**Figure A-12 – NGI Induction Port with Titanium Nitride Coating on Exterior Male Taper**



NGI-0170-6001, Revision B



**Figure A-13 – Three-Piece Pre-Separator (L to R -- Bottom, Insert, Top)**

NGI-0170-6001, Revision B

## **Appendix B**

### **Accessories and Labor Saving Devices for the Next Generation Pharmaceutical Impactor**



**Figure B-1 – Rinsing Caps**

**Left – Male Taper for Induction Port and for Pre-Separator -- 0170-78-0195A-X,**  
**Center -- Face Seal for Straight End of Induction Port-- 4510-78-1042A-C,**  
**and**  
**Right -- Female Taper Plug for Pre-Separator-- 0170-01-0199A-X**